Short Communication Received: June 4, 2012 Revised version: July 12, 2012 Accepted: July 13, 2012

## A Simple HPLC Method for the Determination of Pentyl Gallate and Identification of an Alkaline Degradant

Regina G. KELMANN <sup>1</sup>, Daiane O. DIAS <sup>1</sup>, Ricardo J. NUNES <sup>2</sup>, Helder F. TEIXEIRA <sup>1</sup>, Valquiria L. BASSANI <sup>1</sup>, Cláudia M.O. SIMÕES <sup>3</sup> & Letícia S. KOESTER <sup>1\*</sup>

 <sup>1</sup> Universidade Federal do Rio Grande do Sul, Faculdade de Farmácia, Programa de Pós-Graduação em Ciências Farmacêuticas, Porto Alegre, RS, Brazil.
<sup>2</sup> Universidade Federal de Santa Catarina, Laboratório Relação Estrutura Atividade, Departamento de Química, Florianópolis, SC, Brazil.

> <sup>3</sup> Universidade Federal de Santa Catarina, Laboratório de Virologia Aplicada, Depto. de Ciências Farmacêuticas, Florianópolis, SC, Brazil.

SUMMARY. Pentyl gallate (PG) is a gallic acid n-alkyl ester with a promising antiherpetic topical activity. No method for PG quantification is described in the literature. Therefore, the objective of this investigation is to develop and validate a simple method for determination of PG by High performance liquid chromatography (HPLC) considering its future application in therapy. A nanoemulsion was proposed as PG delivery system due to its low aqueous solubility that limits its delivery through the skin. Chromatographic separation was achieved on a C18 column with 50:49:1 (%, v/v) methanol:water:trifluoroacetic acid (TFA) as mobile phase at a flow rate of 1.2 mL/min; detection at 275 nm. The method was validated for linearity, precision, accuracy, robustness and specificity. The degradant detected under alkaline condition was structurally characterized by time-of-flight mass spectrometry (TOF-MS). All results indicate the method is adequate for PG analysis.

KEY WORDS: HPLC, Nanoemulsion, Pentyl gallate.

918 ISSN 0326-2383

<sup>\*</sup> Author to whom correspondence should be addressed. E-mail: leticia.koester@ufrgs.br