Received: May 11, 2009 Accepted: August 26, 2009

## Technological Development of Hard Capsules of Sertraline Hydrochloride

Adalberto IZQUIERDO C.<sup>1</sup>, Mirna FERNÁNDEZ C.<sup>2\*</sup> & Eddy CASTELLANOS G.<sup>2</sup>

<sup>1</sup> Novatec Laboratories. Ave 25, # 216A03 be/ 216A and 222. Atabey, Playa. Havana City, Cuba. <sup>2</sup> Institute of Pharmacy and Food, University of Havana. Street 23 No. 21425 be/ 214 and 222, La Coronela, La Lisa. Havana City, Cuba.

SUMMARY. Fast release capsules, containing sertraline hydrochloride, pregelatinized maize starch and microcrystalline cellulose were formulated. For this purpose, different technological assays were elaborated being the formulation 2 selected as the better technological variant. Dry powders were filled into hard gelatin capsules. From this formulation were derived capsules with satisfactory technological properties. The quantification of sertraline through HPLC UV detection method was validated for accuracy, precision, linearity and selectivity. The method was linear over the concentration range 0.5 to 0.75 mg/mL and was shown to be highly reproducible. It could be used, without any interference of capsule excipients, for determination of sertraline from solid dosage form. Hard capsules showed an adequate stability during 24 months demonstrating the feasibility from the process of production of this formulation. Parameters, *f1* and *f2*, were used to confirm similarity of dissolution, in deaerated distilled water, of test formulation and capsules of Prosertin as reference product.

KEY WORDS: hard gelatin capsule, sertraline, stability, validation of analytical methods

\* Author to whom correspondence should be addressed. E-mail: mirnafc@uh.cu; mirnafc@yahoo.com