



## Determination of Tolbutamide in Rabbit Plasma by LC–MS/MS and Its Application to Pharmacokinetic Study

Xuebao WANG<sup>1</sup>, Guanyang LIN<sup>2</sup>, Lufeng HU<sup>2</sup>, Xuegu XU<sup>3</sup>,  
Haiyan JIANG<sup>3</sup>, Yibin PAN<sup>3</sup> & Xianqin WANG<sup>1\*</sup>

<sup>1</sup> *Analytical and Testing Center, Wenzhou Medical College, Wenzhou 325035, China*

<sup>2</sup> *The First Affiliated Hospital of Wenzhou Medical College, Wenzhou 325000, China*

<sup>3</sup> *School of Pharmacy, Wenzhou Medical College, Wenzhou 325035, China*

---

**SUMMARY.** A sensitive and selective liquid chromatography–tandem mass spectrometry method (LC–MS/MS) for the determination of tolbutamide in rabbit plasma was developed and validated over the concentration range of 4–1000 ng mL<sup>-1</sup>. After addition of bupivacaine as internal standard (IS), a simplified protein precipitation with acetonitrile was employed for the sample preparation. Chromatographic separation was performed by an Agilent Zorbax SB-C18 column (150 mm×2.1 mm, 3.5 μm). The mobile phase was acetonitrile–1% formic acid in water (50:50 v/v) delivered at a flow rate of 0.3 mL min<sup>-1</sup>. The MS data acquisition was accomplished by multiple reactions monitoring (MRM) mode with positive electrospray ionization (ESI) interface. The lower limit of quantification (LLOQ) was 4 ng mL<sup>-1</sup>. For inter-day and intra-day tests, the precision (RSD) for the entire validation was less than 10%, and the accuracy was within the 94.7% to 105.6% range. The validated method is successfully used to analyze the drug in samples of rabbit plasma for pharmacokinetic study.

---

**KEY WORDS:** LC–MS/MS, Plasma, Pharmacokinetics, Tolbutamide.

\* Author to whom correspondence should be addressed. *Email:* lankywang@yahoo.cn